

Attenuation length measurements of a liquid scintillator with LabVIEW and reliability evaluation of the device^{*}

GAO Long(高龙)^{1,2;1)} YU Bo-Xiang(俞伯祥)^{1;2)} DING Ya-Yun(丁雅韵)¹⁾ ZHOU Li(周莉)¹⁾
 WEN Liang-Jian(温良剑)¹⁾ XIE Yu-Guang(谢宇广)¹⁾ WANG Zhi-Gang(王志刚)¹⁾ CAI Xiao(蔡啸)¹⁾
 SUN Xi-Lei(孙希磊)¹⁾ FANG Jian(方建)¹⁾ XUE Zhen(薛镇)¹⁾ ZHANG Ai-Wu(章爱武)^{1,2)}
 LÜ Qi-Wen(吕绮雯)^{1,3)} SUN Li-Jun(孙丽君)¹⁾ GE Yong-Shuai(葛永帅)^{1,2)} LIU Ying-Biao(刘颖彪)^{1,2)}
 NIU Shun-Li(牛顺利)^{1,2)} HU Tao(胡涛)¹⁾ CAO Jun(曹俊)¹⁾ LÜ Jun-Guang(吕军光)¹⁾

¹⁾ State Key Laboratory of Particle Detection and Electronics, (Institute of High Energy Physics, Chinese Academy of Sciences) Beijing 100049, China

²⁾ Graduate University of Chinese Academy of Sciences, Beijing 100049, China

³⁾ Shanxi University, Taiyuan 030006, China

Abstract: An attenuation length measurement device was constructed using an oscilloscope and LabVIEW for signal acquisition and processing. The performance of the device has been tested in a variety of ways. The test results show that the set-up has a good stability and high precision (sigma/mean reached 0.4 percent). Besides, the accuracy of the measurement system will decrease by about 17 percent if a filter is used. The attenuation length of a gadolinium-loaded liquid scintillator (Gd-LS) was measured as 15.10 ± 0.35 m where Gd-LS was heavily used in the Daya Bay Neutrino Experiment. In addition, one method based on the Beer-Lambert law was proposed to investigate the reliability of the measurement device, the R-square reached 0.9995. Moreover, three purification methods for Linear Alkyl Benzene (LAB) production were compared in the experiment.

Key words: attenuation length, liquid scintillator (LS), Beer-Lambert law, LabVIEW

PACS: 78.20.Ci **DOI:** 10.1088/1674-1137/37/7/076001

1 Introduction

A liquid scintillator (LS) is applied in many neutrino experiments as the target matter. Linear Alkyl Benzene (LAB) is the most important constituent of an LS. Both the Palo Verde neutrino experiment and the CHOOZ neutrino experiment use LAB as the flashing solvent of LS [1]. In the Daya Bay Neutrino Experiment, the new recipe and manufacturing process were adopted to improve the performance of the gadolinium-loaded liquid scintillator (Gd-LS) [2]. It turns out that the performance of the Gd-LS can meet the requirement of the Daya Bay Neutrino Experiment for θ_{13} measurement [3]. In this work, a set-up based on an oscilloscope and LabVIEW was built following an idea that makes it easy to construct as shown in Fig. 1. The LabVIEW program can automatically eliminate the background from every waveform. Otherwise, the accuracy improved due to the use of a filter, the spectral width of light will decrease after it passes through the filter. In the experiment, a

high-precision oscilloscope, a stable high-voltage power supply (HVPS), and an impulse generator were used to improve the precision and stability of the system. A method based on the Beer-Lambert law was proposed to study the reliability of the experimental set-up because there are no samples with an explicit value of attenuation length that can be used to test the device. After comparison with several samples, we took $K_2Cr_2O_7$ solution as the standard sample with KOH (0.05mol/L) as the solvent.

2 The experimental set-up

Besides electronic equipment, the other parts of the experimental facility were placed in a dark room. A light-emitting diode (LED), with peak wavelength of about 430 nm, was used as the light source. The LED was powered by a pulse generator and the pulse width was set to be 4.6 μ s. The light became pointolite after it passed through the collimator, then, it was focused on the fil-

Received 26 September 2012, Revised 15 October 2012

^{*} Supported by National Natural Science Foundation of China (11105160, 11005117)

1) E-mail: gaol@ihep.ac.cn

2) E-mail: yubx@ihep.ac.cn (corresponding author)

©2013 Chinese Physical Society and the Institute of High Energy Physics of the Chinese Academy of Sciences and the Institute of Modern Physics of the Chinese Academy of Sciences and IOP Publishing Ltd

ter and the diaphragm. The aperture diaphragm under the filter was used to get a parallel and very fine beam. The diameter of the light coming out of the aperture was about 2 mm, and it did not change evidently after passing through the sample. The high-voltage power on the photomultiplier tube (PMT:XP2020) [4] was set to be 1500 V. The sample container was a one meter-long cylindrical stainless steel pipe whose capacity was one liter.




Fig. 1. Experimental set up.

3 The experiment

In physics, the attenuation length or absorption length is defined as the distance (λ) in a material where the probability of a particle which has not been absorbed has dropped to $1/e$. With regard to a beam of particles, the attenuation length is the distance in a material when the proportion of incident particles that have not been absorbed has dropped to $1/e$, or when about 63 percent of the particles have disappeared. The particles are

photons produced by the LED in this experiment. The relation between I and I_0 can be expressed as Eq. (1).

$$I = I_0 e^{-\frac{x}{\lambda}}, \quad (1)$$

where I is the light intensity after it passes through the sample, I_0 is the original light intensity, x is the path length of light in the sample, and λ is the attenuation length. The LabVIEW program can read and process about 300 waveforms from the oscilloscope per minute. A waveform consists of ten thousand points, it is transmitted to a computer with a reticle and then handled by the LabVIEW program. Fig. 2 is the data flow diagram of the program to acquire and process the signal. The integral of the waveform is the light intensity. The program takes both ends of the waveform as background dynamically. The statistic of each measurement to a light intensity is more than 3000 and their randomness fits with a Gaussian distribution. For each sample, the intensity of the light changes with the different depths of the sample in the sample container, then the attenuation length can be obtained by fitting the experimental data.

3.1 Stability and precision

The stability of the device which includes the LED, the pulse generator, photomultiplier tube (PMT), HVPS and the data processing system has been studied. As shown in Fig. 3(a), the change of signal is less than 2.5% in 48 hours. The time duration of measuring one sample is about 5 h, so the change is less than 0.3 percent. Therefore, the alterations of the system within an experimental period have little influence on the measurement result. The downswing in Fig. 3(a) can be understood as the PMT will experience tiny and recoverable aging over time. More importantly, as shown in Fig. 3(b), the sigma/mean can reach 0.4 percent after the device is properly adjusted.

3.2 The influence of the filter

The LED is not a monochromatic light source but emits over a significant spectral region. The spectral distribution of the light source used is shown in Fig. 4.




Fig. 2. The data flow diagram of the LabVIEW program.




Fig. 3. The system stability test and data fitting.




Fig. 4. The comparison of the LED's spectra.

The figure also shows that the spectrum width of the beam reduced after it passed through the filter. The following controlled trial will be performed to study how much of an effect the filter will have on the measurements.

In the experiment, we tested one kind of LAB twice. In the first measurement, the attenuation length of the

LAB was 13.96 ± 0.17 m, and the value was 13.93 ± 0.25 m in the second time as shown in Fig. 5. This indicates the measurement system is stable. When the filter was not used, the attenuation length of the LAB was changed to 17.48 ± 0.42 m. It shows that the spectrum's width of test beam will influence the measurement result and the result will become bigger if the spectrum's width becomes wider. The attenuation length of the sample changes along with the wavelength of the light source according to its definition [5]. So the use of a filter makes the results more accurate by narrowing down the spectrum of the light beam. Moreover, the stability of the impulse generator and HVPS is also very important for good results.




Fig. 5. Repeated measurements of LAB and the effect of a filter.

3.3 The reliability evaluation of the device

In this work, a method based on the Beer-Lambert law was proposed to study the reliability of the experimental set-up. The Beer-Lambert law [6] describes the linear relationship between absorbance and concentration of a light-absorbing substance. The general Beer-Lambert law is usually written as:

$$A = -\lg(I/I_0) = \epsilon x c, \tag{2}$$

where A is the measured absorbance, ϵ is the wavelength-dependent absorptivity coefficient, I is the light intensity after the beam passes through the sample, I_0 is the initial light intensity, x is the path length of light in the sample, and c is the concentration of analyte. There are several conditions that need to be fulfilled in order to use the Beer-Lambert law. For instance: making sure the concentration is lower than 0.01M; ensuring that the incident beams are preferably parallel and monochromatic; balancing in chemical equilibria and stabilization in the concentration. If any of these conditions is not fulfilled, deviations will occur from the Beer-Lambert law [7]. The uncertainty of the measured result can be known by measuring the attenuation length of standard samples and

applying the Beer-Lambert law to the measurement results. The standard samples are $K_2Cr_2O_7$ solution which takes KOH (0.05 mol/L) as the solvent. Eq. (3) can be obtained by combining the Beer-Lambert law with the definition of attenuation length where $1/\lambda$ is proportional to the concentration of the standard sample.

$$\frac{1}{\lambda} = \varepsilon \ln(10)c. \quad (3)$$

There is no low limit in concentration, but at very low concentration, the result of absorbance can be erroneous due to the limited resolution of the measuring system or because the signal-to-noise ratio of the light intensity measurement is too low (due to detector noise, photon noise, or the light source fluctuation). Eq. (3) can be applying to the experimental data with a linear fit as shown in Fig. 6. The coefficient of determination (R -square) [8] can be used as one measure of the quality of the linear regression model, it is 0.9995, which means that the measurement device is reliable [9].




Fig. 6. Applying the Beer-Lambert law to standard samples.

3.4 Attenuation length measurements of the Gd-LS

Three different samples including LAB, LS and Gd-LS have been measured in this experiment. They are all used in the Daya Bay Neutrino Experiment. Besides LAB, carboxylic acid, gadolinium chloride, PPO and bis-MSB are used in producing Gd-LS. The Gd-LS prepared using such a recipe has long attenuation length, high photoyield, and long-term stability.

There was a ruler near the tubule which was connected to the sample container, so the position of the sample would be known. The intensity of the outgoing beam was provided by the LabVIEW program. Choosing a position as the reference, we can get the attenuation length with Eq. (1). As shown in Fig. 7, the attenuation length of the LAB is the longest, and the attenuation length of the Gd-LS is 15.10 m.




Fig. 7. The attenuation length measurements of samples.

3.5 Research on purification methods of LAB

In the Phase II project of the Daya Bay Neutrino Experiment, the detector will become bigger, the PMTs will also be placed near the outer wall of each detector to record the light pulses produced when neutrinos pass through the Gd-LS. If the attenuation length of the Gd-LS is too short, the light cannot come out or the intensity of the light will be too weak to be detected by the PMTs. LAB is the main material to synthesize liquid scintillation. The commercially purchased LAB has been purified by vacuum distillation, an aluminum oxide (Al_2O_3) column and a silica gel column, respectively. As shown in Fig. 8, LAB purified by an aluminum oxide column has the longest attenuation length.




Fig. 8. The comparison between different methods in purifying LAB.

Table 1 shows the different attenuation lengths of several samples measured in this work. For $K_2Cr_2O_7$ solution, the unit of concentration is weight percentage.

Table 1. The attenuation lengths of samples.

sample	attenuation length (λ)/m
LAB	19.44±0.62
LS	15.99±0.43
Gd-LS	15.10±0.35
Alcohol	16.82±0.13
0.0365%K ₂ Cr ₂ O ₇	32.42±1.79
0.0611%K ₂ Cr ₂ O ₇	13.04±0.22
0.304%K ₂ Cr ₂ O ₇	3.79±0.04
1.52%K ₂ Cr ₂ O ₇	0.89±0.01
3.04%K ₂ Cr ₂ O ₇	0.46±0.01

4 Conclusions

In this work, the measuring equipment was contrast with LED, filter, PMT, oscilloscope, computer, etc. The LabVIEW program was developed to obtain and han-

dle signal, which made the precision of the whole system reach 0.4 percent. The filter can narrow down the spectrum's test beam width, which enables the test beam to better meet the applicable condition of the formula of attenuation length, therefore the accuracy during measurement will be increased by about 17 percent compared with that without a filter. Besides, the reliability of the measurement system has been evaluated with a method which is based on the Beer-Lambert law. The result shows that the device is reliable. The Gd-LS and the other several samples have been accurately measured in their attenuation lengths. The attenuation length of Gd-LS is 15.10±0.35 m, a Monte Carlo parameter of the Daya Bay Neutrino Experiment for θ_{13} measurement. Moreover, several methods used in purifying LAB have been studied and compared, and it shows that Al₂O₃ is the best candidate.

References

- DING Ya-Yun, ZHANG Zhi-Yong, LIU Jin-Chang et al. Nucl. Instrum. Methods A, 2008, **584**: 238
- DING Ya-Yun, ZHANG Zhi-Yong, ZHOU Peng-Ju et al. Journal of Rare Earths, 2007, **25**(1): 310
- WANG Yi-Fang. Physics, 2007, **36**(3): 207
- Moszynski M. Nucl. Instrum. Methods A, 1993, **324**: 269
- Capone A, Digaetano T, Grimaldi A et al. Nucl. Instrum. Methods A, 2002, **487**: 423
- Baldini F, Giannetti A. Optical Sensing and Spectroscopy, 2005, **5826**: 485
- YIN Hai-Qing. Education Scienc, 2004, **5**(5): 63
- Srivastava A K, Srivastava V K, Ullah A. Econometric Reviews, 1995, **14**(2): 229
- Mahajan M P, Sawant S D. Chen. Tech. Res., 2012, **4**(1): 403